

Heated V-Top Surgery Table

The Heated V-Top Surgery Table is available in two different options, hydraulic base or electric base. This table's integrated patient warming helps prevent hypothermia during small animal surgeries.

Features:

- » controlled heating in temperature ranges from 70°F to 97°F with even heat distribution
- » a two-panel system that adjusts from horizontal to near vertical to hold the patient in the desired position during surgery
- » overall size of the top: 19"W x 58"L
- » each panel contains an 8" x 48" heated area, with an overall heated area of 16" x 48"
- » rotation 360° – electric or hydraulic version
- » the thermostatic control is located on the frame for easy monitoring of temperature
- » two sliding, quick-release, cam cleat tie-down attachments are provided on each side of the attachment rail
- » a 3" stainless steel trough permits fluid runoff into a drain outlet
- » table can be lifted from 30½" to 41" via your choice of electric or hydraulic mechanism
- » weight capacity: 300 lbs for electric and 400 lbs for hydraulic
- » electrical heating system is CSA approved
- » optional table caddy: table base mobility with 2 wheels, 2 wheels are elevated

